



**As every health care worker knows, the risk of needle-stick injuries is a major cause for concern. The Sharp<sub>x</sub>® needle destruction unit has been designed to quickly and easily destroy contaminated needles at the point of use, thereby reducing the risk of needle-sticks without changing accepted procedures or organizational culture.**

### **About the Sharp<sub>x</sub>:**

- **Destroys 18- to 30-gauge hypodermic needles 1/2 to 2 inches in length**
- **One charge destroys as many as 25-100 needles**
- **Needles are destroyed in seconds**
- **Lightweight (2 lbs.) and only 8 1/4" long x 4" wide**
- **Compact charger and NDU fit almost anywhere a 110v wall outlet is available**
- **Built-in drawer opens easily and uses convenient liners for fast debris removal**
- **Efficient handling, easy to transport and use**
- **Unique portal design guides the needle, makes destroying the needle simple**
- **Two visible and one audible alarm keep you informed**
- **Fully recharges in 45-90 minutes**
- **Small size and low cost make it practical**

**Class III FDA Approved (Reg. No. 980015)**

**United States and International Patent Registrations**